

Alternative Methods to Detect COVID-19


A Digital Signal and Image Management
Project

Davoli Sofia - 813479

Filosa Alberto - 815589




PIPELINE



¹ [COVID-19 Cough Recordings | Kaggle](#)

² [Deep Learning and Transfer Learning on COVID-19 | Kaggle](#)



01


AUDIO

Audio Recognition






COVID



Not COVID




AUDIO AUGMENTATION

Only COVID-19 audio



ADD NOISE

Simply add some
random value into
audio



CHANGE SPEED

Change speed of
audio (randomly)




PITCH SHIFTING

Adjust the intonation
of the voice
(randomly)

FEATURE EXTRACTION

- Chroma STFT 1;
- Chroma STFT 6;
- Chroma STFT 12;
- Spectral Centroid;
- Spectral Bandwidth;
- Rolloff;
- Zero Crossing Rate;
- MFCC from 1 to 20.



MODELS



SVM

C: 10

Gamma: 0.0001



LOGISTIC REGRESSION

Solver: LibLinear

Balanced Weights



NEURAL NETWORK

6 Layers

Activation Function: ReLu



CNN


2 Models




BEST MODEL - LOGISTIC

Feature Importance:

- $Coef_{MFCC_3} = -0.57758$
- $Coef_{MFCC_7} = -0.32325$
- $Coef_{MFCC_8} = 0.36815$
- $Coef_{MFCC_{10}} = 0.40355$
- $Coef_{MFCC_{12}} = 0.48419$
- $Coef_{MFCC_{14}} = 0.45921$




	Precision	Recall	F1-Score
Negative	0.92	0.96	0.94
Positive	0.98	0.95	0.96




BEST MODEL (NN)

- Sequential()
- Input(27, 1)
- Convolutional 1D (512, 10)
- Convolutional 1D (256, 10)
- MaxPooling()
- Dense(10)
- Flatten()
- Dense(2) with Softmax



	Precision	Recall	F1-Score
Negative	0.96	0.96	0.96
Positive	0.98	0.97	0.98



LIVE TUTORIAL

You can record your cough to see if you are Positive to COVID-19.

BUT remember, the Accuracy is 97% and it is not a swab!

```
[ ] signal, rate = librosa.load("audio.wav")
    ipd.Audio(signal, rate = rate)
```

/usr/local/lib/python3.6/dist-packages/librosa/core/audio.py:162: UserWarning: PySoundFile failed. Trying audioread instead.
warnings.warn("PySoundFile failed. Trying audioread instead.")

▶ 0:00 / 0:02

```
def kwow_if_you_are_posivite(audio):
    signal, rate = librosa.load(audio)

    #-- Create DataFrame
    new_audio = pd.DataFrame()
    new_audio["filename"] = 0
    new_audio.loc[0] = audio

    #-- Add Features
    new_audio["chroma_stft_1"] = np.mean(librosa.feature.chroma_stft(y = signal, sr = rate)[0])
    new_audio["chroma_stft_6"] = np.mean(librosa.feature.chroma_stft(y = signal, sr = rate)[5])
    new_audio["chroma_stft_12"] = np.mean(librosa.feature.chroma_stft(y = signal, sr = rate)[11])
    new_audio["spectral_centroid"] = np.mean(librosa.feature.spectral_centroid(y = signal, sr = rate))
    new_audio["spectral_bandwidth"] = np.mean(librosa.feature.spectral_bandwidth(y = signal, sr = rate))
    new_audio["rolloff"] = np.mean(librosa.feature.spectral_rolloff(y = signal, sr = rate))
    new_audio["zcr"] = np.mean(librosa.feature.zero_crossing_rate(y = signal))

    for i in range(0, 20):
        new_audio["mfcc_" + str(i)] = np.mean(librosa.feature.mfcc(y = signal, sr = rate)[i])

    y_pred = model.predict(new_audio[new_audio.columns[1:]], verbose=0)
    y_pred = np.argmax(y_pred, axis = 1)

    if y_pred == 1:
        print("Oh Man! You Are Positive to COVID-19! Stay Safe!")
    else:
        print("Lucky You! You are Negative to COVID-19. Wear your mask and Stay Safe!")


    return y_pred
```

```
[ ] kwow_if_you_are_posivite("audio.wav")
```

02

IMAGE

X-Ray Chest Recognition









COVID

Images







Histograms



Images



Histograms





AUGMENTATION

Zooming

From 0% to 20% of zooming




Shearing

Shear angle in counter-clockwise
direction in degrees


MODEL: ResNet - 50

Fine Tuning:

- GlobalAveragePooling2D()
- Dense() with 128 Neurons
- Dropout(0.2)
- Dense(1) with Sigmoid Activation




	Precision	Recall	F1-Score
Negative	0.63	0.79	0.70
Positive	0.85	0.72	0.78




MODEL: EfficientNet – B7

Fine Tuning:

- GlobalAveragePooling2D()
- Dense() with 128 Neurons
- Dropout(0.2)
- Dense(1) with Sigmoid Activation



	Precision	Recall	F1-Score
Negative	0.64	0.76	0.70
Positive	0.84	0.75	0.79





0:00/0:48

CREATED USING
POWTOON



THANKS

Do you have any question?

